



2019-nCoV IgG/IgM RAPID TEST (fingertip blood)

Rapid chromatography immunoassay for reliable detection of 2019-nCoV IgG and IgM antibodies in blood sample

designed for COVID **19**

test result within **15** minutes

98 % correlation with PCR tests

Application

- diagnostic medical device in-vitro
- **no need of unpleasant nasopharyngeal swab**
- uses detection of 2019-nCoV IgG/IgM antibodies in blood
- simple and fast design enables effective testing of large amount of patients

Kit includes

- testing cassette
- tube with extraction buffer
- lancet, dropper
- accessories

Specification

- storage at 2 - 30 °C
- sample stability up to 24 hours at 2 -8 °C
- test result within 15 minutes
- > 98% correlation with PCR method
- control zone for internal procedure control

Rel. sensitivity	96,4 96,9 %
	96,3 %
	96,4 %%

MADE IN EU

